TriVest

#acceptnolimits

TriVest

ESSENTIAL KNOWLEDGE

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Preparing

Pressables

Fix and sprue the wax-up according to the instructions of the ceramics manufacturer and keep a distance of min. 3 mm / 1/8 " to ring base and top.

Crown and bridge alloys

Wax the sprued copings on ring base former and keep a distance of min. 3 mm / 1/8 " to ring base and top.

Spray a thin coat of  **Surfactant** and slightly blow dry.

Mixing/Pouring

- TriVest Liquid press/cast is Frost-sensitive! Temperatures below -5°C / 23°F may cause irreversible cristalization!
- Before mixing, rinse out the clean mixing bowl with water and wipe off.
-> Mixing bowls that are not clean or are kept dry, will withdraw moisture from the investment material!
- Mix the liquid and powder using a spatula for 30 seconds (or use automatic premixing). Then mix for 120 seconds at 350 rpm using a vacuum mixer.
- Available working time: approx. 5–6 minutes (20°C / 68°F, 50 % Liquid concentration).
-> Higher room temperatures will reduce the working time!
- Fill crowns carefully with a fine instrument. Fill the mould ring on the vibrator at the lowest vibration level.
-> Do not vibrate any further after filling!
- Alloy casts – If heating is to be carried out without using a metal ring, remove the rubber sleeve used for investing as soon as possible after complete setting of the investment material (at 20°C / 68°F after approx. 15 minutes); metal mould rings cannot be removed. Pressable moulds must be left to rest and set for min. 15 minutes after being filled. Any deformation caused by moving the mould or early removal from the mould, may cause micro cracks, that can lead to the mould to burst during pressing.
- Do not invest under pressure!
- For speed heating, it is mandatory to stay within the **time window of inserting to furnace 20–25 minutes after first contact of powder and liquid** and the insert temperature of (850°C)!
-> **If you missed this time window, traditional heating process is strongly recommended!**
- Preheat the ring for at least 60 minutes.

Mixing ratio

Only use distilled water for diluting!

Liquid-concentration TriVest Liquid press

Pressables

Inlays / partial crowns	45 – 50 %
Crowns, veneers and bridges	50 – 60 %

100 g TriVest : 22 ml Liquid

Ring size	x 100g-pouch / Qty. Liquid
100 g / 3	1 / 22 ml
200 g / 6	2 / 44 ml
9	3 / 66 ml

Liquid-concentration TriVest Liquid cast

Crown and bridge alloys

precious alloys	%
Inlays / partial crowns	45 – 55 %
Crowns, bridges and primary telescopic crowns	50 – 60 %
Sec. telescopic crowns, CS attachments	55 – 65 %

non-precious alloys	%
Crowns, bridges and primary telescopic crowns	70 – 80 %
Secondary telescopic crowns	75 – 85 %

Density	
TriVest Liquid press	1,2 g/cm ³
TriVest Liquid cast	1,3 g/cm ³

%	45 %	50 %	55 %	60 %	65 %	70 %	75 %	80 %	85 %
TriVest Lq. / H ₂ O	10 / 12 ml	11 / 11 ml	12 / 10 ml	13 / 9 ml	14 / 8 ml	15 / 7 ml	16 / 6 ml	17 / 5 ml	18 / 5 ml

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Concentration du TriVest Liquid cast

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Safety instructions

- Investment materials contain quartz. Do not inhale dust! Risk of damage to lungs (silicosis, lung cancer). Recommendation: Use gas mask type FFP 2 – EN 149:2001. Open bag along the perforation or cut open with scissors and avoid generation of dust when filling into mixing bowl. Rinse out empty bag with water before crumpling up. Following eye contact: Remove contact lenses, hold eyelids open and rinse thoroughly with clean, running water for at least 15 minutes. Consult an eye specialist.
- Remove dust at workplace only when moist.
- To avoid dust during deflasking, place the mould in water after it has completely cooled down after casting until it is thoroughly moistened.
- Use suction extraction system with fine dust filter when blasting.
- Remove dried-up mixing liquid only when moist (fine dust). Seal bottle securely after use.
- The setting reaction of the investment-material generates heat. To avoid injuries, please use gloves or suitable muffle tongs.
- Combustion gases resulting from the preheating of casting moulds must be discharged via a suitable smoke outlet.
- Place all moulds in the furnace within 10 seconds – then keep the furnace door closed for 15 minutes!
- Never touch hot moulds with your hands. Use suitable mould tongs to avoid burns.
- After casting, allow the moulds to cool down until warm to the touch in a protected location!

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Contains: Quarz (SiO₂), Cristobalit
H372: Causes damage to lungs through prolonged or repeated exposure: Inhalation
P260: Do not breathe dust/fumes/gas/mist/vapours/spray.
P285: In case of inadequate ventilation wear respiratory protection.
P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

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General instructions

- The ideal processing temperature for TriVest is 20 °C. To keep it constant at a higher ambient temperature, use a temperature control cabinet if necessary. At a high ambient temperature also place the mixing bowl and liquid in the temperature control cabinet.
- Observe shelf life dates of investment material and mixing liquid! Do not use without checking after the specified date.
- Do not use crystallized mixing liquid any longer. The higher the concentration of mixing liquid, the greater the expansion.
- Data on the concentration of the mixing liquid are based on our own experience and tests and are to be regarded as a reference value only.

Application-related recommendations made by us, whether given verbally, in writing or by way of practical instructions, are based on our own experience and tests and may therefore be regarded only as general guidelines.

Our products are subject to continuous development. We therefore reserve the right to change the design and composition of our products. Work results may be influenced by such parameters as temperature, mixing intensity, investment system, etc.

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HPdent GmbH