Print-Parameter InovaPrint

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Print-Parameters InovaPrint

PRINTING PARAMETERS RECOMMENDATION

Depending on the type of printer and individual expertise, the parameters can be individually adjusted and optimised. Our **InovaPrint** resins are compatible with all open source LCD and DLP in the range of 385 – 420 nm.

PREPARATIONS BEFORE STARTING TO PRINT

1. **InovaPrint** resins should be used at a minimum room temperature of 21 °C (70°F).
2. Shake the bottle well for 2 min. before each use.
3. Important: observe the general precautions when using chemicals! Always use nitrile gloves and safety glasses.
4. After shaking the bottle, leave the resin for 10 min. to let air bubbles get out of the resin.
5. Excess resin can be poured back into the bottle using a special 3D printing resin filter.
6. Prevent dust or other dirt particles from entering the resin VAT.

GENERAL SUGGESTIONS

for getting a good printing result:

- ▶ The more volume (weight) the print object has, the more burn in layers are suggested to be used.
- ▶ We suggest to hollow your print objects always – wall thickness of 2,5 mm is suggested to decrease the adhesion with the FEP film.
- ▶ Always keep in mind:
The larger the contact area on the FEP film per layer, the greater the risk of failures due to the high adhesion not releasing the object.
The larger the contact area, the larger the „lift-off height after printing in mm“ value must be.
- ▶ We recommend the cleaning process with **InovaPrint wash** to work safely, quickly and economically and to leave no white residues on the objects.
- ▶ Cure your models in a high power curing chamber for approximately 20–30 minutes at 60° C (140°F) – The preferred curing wavelength is between 385–420 nm.
- ▶ If needed, cure **InovaPrint ocean | blue** in a high power curing chamber with max. 80 % power for approximately 10 minutes at 30° C (86°F) – The preferred curing wavelength is between 385 – 420 nm.

Our InovaPrint resins are validated with  ACKURETTA and  ASIGA®.

The print parameters are integrated in the libraries of the slicing softwares.

* Please follow printer manufacturers instructions to actualize your parameters library..

Printer Technology „Standard RGB LCD Panel“

For example: Ackuretta Freshape 120, Dentiq, Anycubic Photon/Phroton S, Phrozen

Colors/Resins		Layer Thick-ness μm	Base Layers Exposure Time sec.	Number of Base Layers	Base Exposure Off Time sec.	Exposure Off Time sec.	Layers Exposure Time sec.	Lift Height mm	Lift Speed mm/min.	Lowering Speed mm/min.
3D PRINTED MODELS										
HPdent recommends a resolution of 50 μm to print a high quality/precise model (print-results will be more precise or at least equivalent to class IV plaster model)										
	InovaPrint orca black opaque model	50	140	8	2	2	14,5	8	130	100
	InovaPrint beluga white opaque model					1	13,5			
	InovaPrint dolphin grey opaque model									
PRINT 'N PRESS OR CAST										
	InovaPrint ocean blue pressable castable	50	130	5	2	1,5	10,5	6	130	100

Printer Technology „Monochrome LCD Panel“

For example: Anycubic Photon Mono, Phrozen Mono etc.

Colors/Resins		Layer Thick-ness μm	Base Layers Exposure Time sec.	Number of Base Layers	Base Exposure Off Time sec.	Exposure Off Time sec.	Layers Exposure Time sec.	Lift Height mm	Lift Speed mm/min.	Lowering Speed mm/min.
3D PRINTED MODELS										
HPdent recommends a resolution of 50 μm to print a high quality/precise model (print-results will be more precise or at least equivalent to class IV plaster model)										
	InovaPrint orca black opaque model	50	20 – 50	8	3	1 – 2	2	8	60	60
	InovaPrint beluga white opaque model									
	InovaPrint dolphin grey opaque model									
PRINT 'N PRESS ODER CAST										
	InovaPrint ocean blue pressable castable	50	50 – 60	5	3	1	3	6	90	150

Printer Phrozen Sonic Mini 8 K

Farben/Resin	Layer Thick- ness μm	Base Layers Exposure Time sec.	Number of Base Layers	Base Exposure Off Time sec.	Exposure Off Time sec.	Layers Exposure Time sec.	Lift Height mm	Lift Speed mm/min.	Lowering Speed mm/min.	
3D PRINTED MODELS										
HPdent recommends a resolution of 50 μm to print a high quality/precise model (print-results will be more precise or at least equivalent to class IV plaster model)										
	InovaPrint orca black opaque model	50	35	7	4	1–2	5	8	60	150
	InovaPrint beluga white opaque model									
	InovaPrint dolphin grey opaque model									
PRINT 'N PRESS ODER CAST										
	InovaPrint ocean blue pressable castable	50	60	5	3	1	5	6	90	150

Elasticity modulus:	2,1 GPa
Flexural strength:	82 MPa
Tensile modulus:	1,9 GPa
Tensile strenght:	55 MPa
Shore D hardness:	84
Water sorption:	0,38 %

Workflow – 2-stage cleaning with InovaPrint wash

1. Fill both containers (Prewash; Cleanup) with InovaPrint wash.


2. Place the 3D objects in the prewash container and clean them with the ultrasonic cleaner for 2 to 10 minutes, depending on their size.

3. The liquid has removed excess resin. Remove the 3D objects ...

4. ... and place them in the cleanup container. Clean here for another 2 to 10 minutes with the ultrasonic cleaner.


5. After cleanup, remove the 3D objects, rinse with water and light cure.

6. Great results and no white residue on the objects.

7. The used InovaPrint wash, which is now saturated with resin, can easily be recycled with UV light: Simply leave the container in daylight for a day.

8. The resin will settle on the bottom of the container. The clear liquid on top can now be carefully removed with a syringe and used for the next workflow.